Eli-VET™

Designed for veterinary market

Blood Glucose Monitoring System

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Contents

- I. Important Safety Instructions
- II. Operating Instructions Before First Use
- III. Environmental Conditions Required for System Use (Normal Operating Conditions)
- IV. Measurement Procedures for Evaluating System Performance Characteristics
- V. Recommended Measures for Preventing Infection Before Using the device
- VI. Main Structure of the Product
- VII. Sample Types Used for Calibration
- VIII. Description of Symbols Used in Labels and User Guide
- IX. Operating Procedures to Follow When Using the Meter
- X. Description and Interpretation of Symbols Displayed on the Meter Screen
- XI. Reading Meter Memory
- XII. How to Set Time and Voice
- XIII. Identification of Error Messages During Meter Operation and Recommended Measures
- XIV. Types of Samples Used and Collection Methods
- XV. Recommendations for Actions Based on Measurement Results
- XVI. Specifications of the Product
- XVII. Scope of Application
- XVIII. Contraindications
- XIX. Principles of the Testing Method
- XX. Factors Affecting the Test and Precautions
- XXI. Production date and Shelf Life
- XXII. List of Accessories
- XXIII. Maintenance, Quality Control Testing, and After-sales Service of the Meter
- XXIV. Transport and Storage of the Device
- XXV. Electromagnetic Compatibility Description
- XXVI. Reference Range

Customers should read all the information in the user guide before testing. All illustrations and texts in this guide are for illustrative purposes only; refer to the device itself for specifics.

I. Important Safety Instructions

DANGER - Improper use of electrical device can cause electrocution, burns, fire, and other hazards. Always perform basic safety precautions, including those listed below. Read the followings before using the device.

1. Do not immerse the device in liquid or place it where it could fall into liquid.
2. Use the device only as described in the user guide.
3. Do not use accessories not provided or recommended by the manufacturer.
4. Do not use the device if it is not working properly or if it has been damaged.
5. Do not allow the device or its cord to come into contact with surfaces that are too hot to touch.
6. Do not place any items on top of the device.
7. Unless specifically required by the user guide, do not allow any objects to fall into or be inserted into any openings, ducts, or seams of the device. Do not use the device where aerosol spray is being used or where oxygen is being administered.
8. If the device is out of use due to maintenance or other treatments, disinfect it before using it again.
9. Keep these instructions

II. Operating Instructions Before First Use

1. Battery Installation and Usage Instructions

This product requires two AAA batteries. When the "battery" icon appears, it indicates that the battery power is low, and the batteries should be replaced promptly. The battery installation steps for the device are as follows:

- a) Before installing the batteries, check if there are usable batteries in the device and use the meter normally as indicated in the user guide if there are.
- b) If there are no batteries in the meter or if the battery power is insufficient, flip the meter over with the back facing up and open the battery cover by pressing upward with your finger (Note: when the battery power is insufficient, remove the batteries that are low on power).
- c) Install the new batteries correctly according to the polarity indications.
- d) Close the battery cover and ensure it is securely fastened.

Note: If the meter is not going to be used for a long time, please remove the batteries. For safety, keep the meter and its accessories out of reach of children.

Dispose of used batteries, meter, and accessories properly as they can pollute the environment. Do not dispose of them indiscriminately. Follow local regulations for disposing of used batteries, meters, and accessories.

2. Instructions for Connecting Accessories or Other Devices

When using this meter for blood glucose testing, it is necessary to use it in conjunction with Eli-VET™ test strips and Eli-VET™ control solution. Additionally, use it in conjunction with a lancet and lancing device. Refer to Eli-VET™ test strip user guide for detailed usage.

III. Environmental Conditions Required for System Use (Normal Operating Conditions)

1. Environmental Temperature: +5°C(41°F) to +40°C(104°F);
2. Relative Humidity: ≤80%;

3. Atmospheric Pressure: 860hPa to 1060hPa;
4. Device Power Source: DC 3.0V (Two 7# AAA batteries).

IV. Measurement Procedures for Evaluating System Performance Characteristics

Evaluate the performance of the testing system using Eli-VET™ control solution, and ensure that the test results fall within the control range of the control solution.

V. Recommended Measures for Preventing Infection Before Using the Meter

1. Ensure that the lancet is being used for the first time.
2. If you come into contact with pet blood, use protective gloves or other protective measures.
3. If you have doubts about the operation and need training, please contact the manufacturer.


VI. Main Structure of the Product

1. Structural Composition


The meter consists of a detection module, data processing module, display module, signal output section, power circuit, and key control circuit etc.

Note: When using this product for blood glucose testing, it must be used in conjunction with Eli-VET™ test strips. Additionally, a lancet and lancing device are also needed for use.

2. Names and Descriptions of Components

Front of the Meter


Back of the Meter

Main Button Functions: Power on, power off, assist in selecting animal type.

Settings Button 'S' Functions: Time and voice settings; assists in reading blood glucose memory values.

Memory Button 'M' Functions: Reads blood glucose memory values; assists in setting time and voice.

3. Description of the Meter Display Screen

4. Blood Glucose Test Strip Instructions

Blood Glucose Test Strips: Eli-VET™ Test Strips (Glucose Dehydrogenase Method):

VII. Types of Samples Used for Calibration

Fresh capillary whole blood or venous whole blood.

VIII. Description of Symbols Used in Labels and User Guide

1. "—" Direct current symbol.
2. "⚠" Caution, danger. Consult the document in all cases where this symbol is marked.
3. "☠" Biohazard warning.
4. "+" represents the positive terminal of the battery, "-" represents the negative terminal.

IX. Operating Procedures to Follow When Using the Meter

(1) Measure blood glucose using the normal operating mode of the meter.

(Before using, please confirm that you are using Eli-VET™ Test Strips, avoid using the wrong strips.)

(All illustrations are for example purposes only, refer to the actual product.)

Step 1: Read the Calibration Card

Place the calibration card in the corresponding area of the meter.

(Note: Usually, one calibration card is only compatible with one box of test strips shared by cats and dogs.

Each time you use the next box of test strips, the meter should be recalibrated with the new calibration card from the new packaging.)

Step 2: Insert the Test Strip

1. Open the test strip packaging and take out a blood glucose test strip. (Note: Once the sealed bottle is first opened, store as instructed, and it remains valid for 90 days. Immediately cap the bottle after removing a test strip to prevent the remaining strips from getting damp.)

2. Hold the middle part of the test strip (Note: Do not touch the reaction area of the strip), and align it with the strip port.

3. Electrically insert the blood glucose test paper into the blood glucose test paper port (the test paper can be inserted in the off and on state). The meter will display a calibration code and then voice prompt for blood collection.

(Note: When the meter screen displays the corresponding blood drop symbol, it indicates that the strip is properly inserted. Ensure the front of the test strip faces the same direction as the meter screen.)

When using the meter, please confirm that the calibration code displayed on the glucometer matches the calibration code on the test strip packaging.

If you are not clear about the calibration code, you can pull out the test strip and re-insert it. Inconsistency in the calibration code can lead to significant measurement errors.

(Note: The calibration codes shown in the diagrams are for illustrative purposes only, use the actual calibration codes for dogs and cats provided with the test strips.)

Step 3: Select Animal Type

1. After inserting the test strip, the last default test animal type and corresponding calibration code will be displayed. The meter defaults to dog mode.

2. If the current test subject is a cat or dog, press the "Main Button" to switch the animal type selection.

(Note: The currently selected animal type will be remembered for the next test startup, reducing unnecessary

repetitive selections.)

Step 4: Blood Collection


Note: Blood collection can be done using a disposable lancet with a safety lock or a lancing device with a compatible lancet. A disposable lancet with a safety lock does not require a lancing device. The following blood collection method uses one of these options.

I. Blood collection using a lancing device and compatible lancet

1. Take out the lancing device, unscrew the cap (the first + second sections are the cap part).
2. Take out the lancet, insert the long end of the lancet into the pen core (cylindrical groove).
3. Twist the protective cap of the lancet two to three turns and then pull it off.
4. Place and tighten the cap of the lancing device.
5. Rotate the first and second sections of the lancing device cap to adjust the needle penetration depth. 1→5, from shallow to deep, adjust according to the animal's blood collection site. Generally, set to position 3, thicker sites are recommended to be adjusted to >3. The number pointed to by the arrow represents its depth.
6. Finally, pull the pen plunger (the longer part at the end of the lancing device) until the middle button pops up.

7. Take out the prepared lancing device.

(Note: The lancing device is for individual animal use only and should not be shared. Lancets are for single use only.)

8. Press the tip of the lancing device vertically against the blood collection site as shown in the diagram below. (Tip: If needed, lightly rub the blood collection site or warm it with a warm cloth before collecting blood to improve blood sample acquisition.)

Ear Margin Vein

Elbow (Dog)

Paw Pad

Inner Lip (Dog)

9. Press the button on the lancing device, causing the needle to prick.
10. Then gently press around the wound to form a sufficient amount of blood droplet. (Tip: Avoid contaminating the blood droplet to prevent measurement errors.)

II. Blood Collection Using a Disposable Lancet with a Safety Lock

1. Take out the prepared disposable lancet with a safety lock. (Tip: Lancets are for single use only.)
2. Twist the protective cap, pull it out, and discard.
3. Align the lancet with the blood collection site and press down quickly.
4. Then gently press around the wound to form a sufficient amount of blood droplet.

(Tip: Avoid contaminating the blood droplet to prevent measurement errors.)

Step 5: Absorb Blood

1. When the corresponding blood droplet symbol on the screen flashes, gently touch the blood droplet to the top of the test strip to fill the reaction area with blood in one go.

(Tip: ● Do not drop blood on the surface of the test strip. ● Do not interrupt the absorption process. ● Do not absorb blood repeatedly.)

2. After the blood is fully absorbed, place the meter gently on the table and wait for the measurement result.

(Tip: Do not hold the meter in your hand for a long time to avoid causing an increase in body temperature.)

Step 6: Display Results

Approximately 1uL of blood sample is needed for the blood glucose test. After the countdown, the test result is displayed, measurement unit: mmol/L (mg/dL in USA) .

Step 7: Power Off

Press the button to eject the test strip, and the meter can automatically power off.

(Tip: After testing, collect the used test strips and dispose of them in a "medical waste bin". After each test, properly handle the used lancet and blood glucose test strip. Used test strips and lancets are medical waste, so be sure to comply with local regulations on the disposal of animal medical waste and do not reuse them. After the test is complete, immediately eject the test strip to avoid accelerating battery consumption of the meter.)

(2) Maintenance After Using the Meter

Do not wet the meter. There should be no liquid contamination at the meter test strip port or inside. If there are stains on the surface of the meter, wipe them with a soft cloth. Store in a dry and normal temperature environment. You can use a small amount of 75% alcohol for meter disinfection. Avoid placing in strong electromagnetic environments. Do not drop or forcefully operate the device, as it may cause damage and render the device unusable.

(3) If the meter is not in use for a long time, please remove the batteries.

(4) Data Units Displayed by the Meter

Blood Glucose: mmol/L (millimoles per liter). (mg/dL in USA) (1mmol/L=18mg/dL)

X. Description and Interpretation of Symbols on the Meter Display

1. When the test strip is inserted and the device is powered on, the current calibration code is displayed, for example, the calibration code "588" for a dog. The calibration code displayed on the meter must match the calibration code indicated on the test strip packaging.

2. When the "Test Strip" icon on the display blinks, it waits for the user to insert a test strip.

3. When the water droplet-shaped blood sample symbol on the display blinks, it waits for the user to add a blood sample.

4. After adding a blood sample, the meter enters the timing test mode.

5. The value that appears after the countdown is the result of this test. When testing the control solution, the control solution icon is displayed on the screen.

Blood glucose unit: mmol/L (millimoles per liter) (mg/dL in USA)

6. When the "Battery" icon appears, it indicates that the battery power is insufficient, and the batteries should be replaced promptly.

7. The date and time displayed at the top of the screen record the date and time of this test.

8. To confirm that the display is working correctly, hold down the main button "  " when the device is off, and the full display should be visible. If it is not fully displayed, please contact the after-sales service department of our company.

XI. Reading Meter Memory

With the device off, turn it on by pressing the main button, and the interface displays "OK". When the "OK" interface is displayed, press the memory key "M" to enter the historical test result browsing. Press the memory key "M" to select the item to browse (7-day average value for dogs, 14-day average value for dogs, 28-day average value for dogs, 7-day average value for cats, 14-day average value for cats, 28-day average value for cats, dog glucose value, cat glucose value, control solution test value), and press the setting key "S" to view the result value. The specific operations are as follows:

The first time you press the memory key "M", the screen displays the 7-day average value of dog glucose levels;

The second time you press the memory key "M", the screen displays the 14-day average value of dog glucose levels;

The third time you press the memory key "M", the screen displays the 28-day average value of dog glucose levels;

The fourth time you press the memory key "M", the screen displays the 7-day average value of cat glucose levels;

The fifth time you press the memory key "M", the screen displays the 14-day average value of cat glucose levels;

The sixth time you press the memory key "M", the screen displays the 28-day average value of cat glucose levels;

The seventh time you press the memory key "M", the screen displays the last tested dog glucose value, press the setting key "S" to sequentially view 500 sets of dog glucose result values;

The eighth time you press the memory key "M", the screen displays the last tested cat glucose value, press the setting key "S" to sequentially view 500 sets of cat glucose result values;

The ninth time you press the memory key "M", the screen displays the last tested control solution test value,

press the setting key "S" to sequentially view 10 sets of control solution test values.


The above items cycle through. Long press the memory key "M" until the "OK" character appears to exit the memory value browsing and average value query function.

The top-left corner of the screen shows the test time (month, day, hour, minute) of the current memory value, and the bottom-right corner shows the sequence number of this memory value or the number of days for the average value. When browsing average values, dog, cat, or control solution test values, the corresponding icon is displayed above the result value.

XII. How to Set Time and Voice

With the device off, press the main button to turn it on, and the interface displays "OK". When the "OK" interface is displayed, press the setting key "S" to enter the time and voice settings. Press the setting key "S" to sequentially select the item you want to adjust (minutes, hours, days, months, years, voice), and press the memory key "M" to change the value of that item. The voice can be turned on or off. After completing the settings, long press the setting key "S" until the "OK" character appears, exiting the time and voice setting function.

XIII. Identification and Recommended Actions for Error Messages During Meter Operation

1. "EXX" - This symbol indicates an error has occurred, e.g., E07

Error Indicator	Description	Possible Reasons
E07	Voltage at the test strip end exceeds the range when the device is turned on	Wrong test strip inserted; invalid test strip used
E08	Voltage at the test strip end exceeds	Wrong test strip inserted; invalid test strip

	the range after inserting the test strip when prompted to insert	used
E09	Test strip pulled out during the test process	Operational error
E12	Calibration card not detected or calibration card malfunction	Calibration card not detected or calibration card malfunction
E13	Sample error	Sample error
E14	Impedance value exceeds the range error	Test strip moved during the test process or test strip malfunction
E16	Wrong card detected	Calibration card found incompatible with the meter
E17	Saved parameters in the meter are incorrect	Meter malfunction

Other possible reasons: Incorrect operation method (steps); test strip and meter mismatch; Insufficient blood sample; blood drawn before inserting the test strip; blood not drawn when prompted to draw blood; test strip exposed to air for an extended period; finger touched the reaction area; meter malfunction, etc.

Recommendations: a) Carefully read the instructions and retake the test with a valid test strip using the correct method.

b) Please contact our company's after-sales service department.

2. "Err" - This symbol, along with a "thermometer" icon, indicates that the surrounding environmental temperature is outside the normal operating conditions. The meter should be used in an environment between +5°C (41°F) and +40°C (104°F).


3. "code exp" - This symbol indicates that the calibration card has expired, and the test strips in the corresponding test strip box or test strip bottle are expired. A new calibration card from a new package or bottle should be used.

4. "Lo" - This symbol indicates that the blood glucose measurement is below 1.1mmol/L (20mg/dL) .

Recommendations: a) Check if the measurement operation is correct and retake the test with a valid test strip. If the results are still consistent or approximately the same, visit a veterinary clinic. When symptoms and test results are abnormal, follow the veterinarian's guidance.

b) Please contact our company's after-sales service department.

5. "Hi" - This symbol indicates that the blood glucose measurement is above 33.3mmol/L(600mg/dL) .

Recommendations: a) Check if the measurement operation is correct and retake the test with a valid test strip. If the results are still consistent or approximately the same, visit a veterinary clinic. When symptoms and test results are abnormal, follow the veterinarian's guidance.

b) Please contact our company's after-sales service department.

Please note: Promptly remove the used test strip after the measurement to prevent loss of memory values when turning the device on again.

XIV. Sample Types and Collection Methods

This system tests fresh capillary whole blood or venous whole blood from cats or dogs. Blood collection sites can be selected from the ear edge vessels or the side of the paw for cats or dogs, the elbow callus or inner lip for dogs.

Fresh Capillary Whole Blood Collection Method: Use a safety-lock cap disposable peripheral blood needle or a blood collection pen with a disposable sterilized blood needle to puncture the animal's blood collection site.

Venous Whole Blood Collection Method: Venous whole blood samples can be fresh heparin anticoagulated whole blood, EDTA anticoagulated whole blood, or potassium oxalate anticoagulated whole blood. Samples should be used within 10 minutes of collection by professional veterinary personnel. Hemolyzed or coagulated samples may lead to inaccurate test results.

Note: Only use fresh peripheral blood samples or venous whole blood. Do not use arterial blood, serum, or plasma samples.

XV. Recommendations for Users Based on Measurement Results

1. When the results obtained with this test system do not match the condition of the animal, please check whether the operation method is correct and retake the test with a valid test strip. If the situation remains the same, consult a veterinarian promptly. Do not change medication or medical advice solely based on the test results without consulting a veterinarian.

2. When the measurement results fall outside the analytical range, the system will prompt the user with a symbol display.

Blood Glucose: "Lo" indicates that the blood glucose measurement is below 1.1mmol/L (20mg/dL); "Hi"

indicates that the blood glucose measurement is above 33.3mmol/L (600mg/dL).

3. Regularly compare the blood glucose measurements with those from a pet clinic. When symptoms and test results are abnormal, follow the veterinarian's guidance.

XVI. Product Specifications

1. Meter model specification: Eli-VET™ Glucose meter
2. Meter performance parameters
 - a) Power supply: DC 3.0V (two 7# AAA batteries);
 - b) Display method: Blood glucose testing: mmol/L (mg/dL in USA)) display and test time display;
 - c) Test result memory count: 500 memory entries for dog blood glucose, 500 for cat blood glucose, and 10 for quality control liquid tests;
 - d) Dimensions (mm): (115±5) × (59±5) × (20±5);
 - e) Total weight (g): 65±10 (excluding batteries);
 - f) Maximum power (mW): 800.
3. Specification and model of applicable test strips: Blood glucose test strips: Eli-VET™ Test Strips (Glucose Dehydrogenase Method):

XVII. Application Scope

This product is used in conjunction with Eli-VET™ glucose test strips for the in vitro quantitative measurement of glucose concentration in peripheral and venous whole blood in animals (cats and dogs). It is suitable for pet owners at home or veterinarians in medical institutions to monitor the blood sugar of cats and dogs. It is only for monitoring the effectiveness of blood sugar control in cats or dogs, not for diagnosing or screening diabetes in cats or dogs, nor can it be the basis for adjusting treatment drugs for cats or dogs.

XVIII. Contraindications

No contraindications have been identified.

XIX. Principle of Detection Method

The detection is based on the principle of enzyme electrodes and dry chemistry technology. Glucose in the blood reacts with glucose dehydrogenase fixed on the test strip, causing electron transfer and generating a current. The magnitude of the current is directly proportional to the glucose concentration.

XX. Factors Affecting the Test and Precautions

1. Factors Interfering with Test Results due to Physiological Diseases
Hematocrit is greater than 70%.
2. Factors Affecting Test Results due to Improper Operation
 - a) Insufficient blood volume, not fully filling the reaction area of the test strip.
 - b) Insufficient initial blood sample volume, and then re-adding a blood sample for testing. Ensure that the blood sample volume for each test strip is sufficient, and the same test strip cannot be used to add blood samples twice.
 - c) Blood sample contamination (blood sample obtained by forced squeezing or blood sample containing bubbles).
 - d) Removing the test strip after adding blood and then re-inserting it into the meter.
 - e) Using serum or plasma samples. Only suitable for whole blood testing, do not use serum or plasma samples.

- f) Reusing test strips. Blood glucose test strips are for in vitro testing only and are for single-use.
- g) Testing with test strips stored or refrigerated in an environment below 4°C(39°F) or frozen. Using test strips stored at high temperatures above 35°C(95°F) . Improper storage of test strips will affect the test results. To ensure the best shelf life of the test strips, store them in a dark and airtight environment at 4°C (39°F) to 35°C(95°F) , do not freeze the test strips or store them in places with high temperature and humidity.
- h) Exposing the test strip to the environment for a longer time than required after opening the packaging. The test strip should be used within 3 minutes of being taken out in an environment of 5°C (41°F) to 40°C (104°F) and relative humidity of 10% to 80%. In high humidity environments (relative humidity >80%), the test strip should be used immediately after being taken out.
- i) Using expired blood glucose test strips or calibration cards. Do not use test strips and calibration cards that have expired.
- j) Not allowing the meter and test strips to fully equilibrate with the ambient temperature of the usage environment. If there is a significant difference between the storage temperature of the test strips and the meter and the temperature of the testing environment, allow the test strips and the meter to equilibrate in the testing environment for more than 30 minutes before testing.
- k) Use of incorrect test strip or control solution models, such as those from other parameters, manufacturers, or models. Eli-VET™ pet glucose test strips, pet glucose meters, and pet glucose control solutions should be used as a set. Do not replace them with products from other brands.

3. Note

Thank you for choosing the Eli-VET™ glucose monitoring system. This product is for in vitro quantitative measurement of glucose concentration in peripheral and venous whole blood in animals (cats and dogs). It helps you monitor your pet's blood glucose levels and is suitable for pet owners at home or veterinarians in medical institutions for monitoring the blood glucose of cats or dogs. It is only for monitoring the effectiveness of blood glucose control in cats or dogs, not for diagnosing or screening diabetes in cats or dogs, nor as a basis for adjusting treatment medications. The testing frequency and testing period vary for different animals, and should be conducted in accordance with medical advice. Testing should be done under the guidance of a professional veterinarian, and interpretation of test results should be determined by a professional veterinarian, not by the user. This product has undergone quality control and is approved for distribution only after testing confirms no errors.

Customers should use the device according to the usage described in the guide. If the device is not used according to the manufacturer's specified method, the protection provided by the device may be compromised. If the device is not working properly or has been damaged, do not use it.

Use control solution to perform quality control of the testing system if you suspect the device or test strip is not in a normal condition; if you doubt the test results or if the results are inconsistent with usual conditions; or when the device is damaged.

XXI. Production date and Shelf Life

Production date is indicated on the device label.

Shelf life: 5 years if stored at normal room temperature.

XXII. Accessory List

1. Glucose meter: Eli-VET™ glucose meter x 1
2. User guide x 1
3. Batteries: AAA batteries x 2

Specific accessory lists are subject to actual packaging.

XXIII. Maintenance, Quality Control Testing, and After-sales Service of the Meter

1. Maintenance

Do not get the meter wet, and ensure that there is no liquid contamination at the meter's port or interior. If the surface is stained, clean it with a soft cloth. Wipe and disinfect the meter's surface with 75% alcohol. Store in a dry place at room temperature. Avoid dropping the meter. If the meter is not to be used for a period of time, it is recommended to remove the batteries.

If the device is to be stopped for repair or handling, it should be disinfected first to eliminate or reduce the risks associated with cessation of use, transportation, or handling.

Apart from batteries, other parts of the meter should only be inspected or provided by our company.

2. Quality Control Testing

Quality control testing ensures the meter and test strips are functioning correctly, and operations are accurate, thus guaranteeing the precision of test results. Conduct quality control testing in the following situations:

- a) When you suspect the meter or test strips are not functioning normally.
- b) When you doubt the test results or if they are inconsistent with usual conditions.
- c) When the meter is damaged.

3. Elimedical warrants the product will be free from defects in materials and workmanship for 5 years from the date of original purchase. This warranty is limited to replacement due to defects in parts or workmanship. Elimedical shall not be required to replace any units that malfunction or are damaged due to abuse, accidents, alteration, misuse, neglect, or failure to operate the instrument in accordance with instructions. Elimedical assumes no liability for malfunction or damage caused by the use of control reagents or test strips other than Elimedical ones. In no event shall Elimedical be liable for indirect, special or consequential damages, even if Elimedical has been advised of the possibility of such damages.

XXIV. Transportation and Storage of the Meter

1. Transportation

During transportation, the meter should be kept at a temperature of $-40^{\circ}\text{C}(-40^{\circ}\text{F})$ to $+55^{\circ}\text{C}(131^{\circ}\text{F})$ with a relative humidity not exceeding 93%. It should be properly arranged according to the markings on the packaging box, and protected from heavy pressure, collision, severe vibration, and direct rain or snow. Other transportation requirements should be implemented according to the stipulations in the order contract.

2. Storage

The packaged meter should be stored indoors, where the temperature is between $-40^{\circ}\text{C}(-40^{\circ}\text{F})$ and $+55^{\circ}\text{C}(131^{\circ}\text{F})$, relative humidity does not exceed 93%, there are no corrosive gases, strong mechanical vibrations are absent, and the environment is clean, hygienic, and well-ventilated.

XXV. Electromagnetic Compatibility Information

1. Assessment of Electromagnetic Environment

It is advisable to assess the electromagnetic environment before using the device.

2. Prohibitions

The meter must not be used near strong radiation sources, such as unshielded RF sources, as it may interfere with the normal operation of the meter. In dry environments, especially those containing synthetic materials (synthetic fabrics, carpets, etc.), the use of this device may lead to destructive electrostatic discharges, resulting in erroneous conclusions.

XXVI. Reference Ranges

Dog glucose range: 3.3-6.7mmol/L (59-121mg/dL)

Cat glucose range: 3.5-7.5mmol/L (63-135mg/dL)

In Canada and UK, the standard unit of blood glucose measure is mmol/L. In the U.S., meters report results in milligrams of glucose per deciliter of blood, or mg/dL.

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Importer in USA: ELIMEDICAL-US LLC.