



Small Syringe.
Meaningful Care.



November 2026

Pet Diabetes Awareness Month

Every November, Pet Diabetes Month provides an opportunity to raise awareness about diabetes mellitus, a common endocrine disorder affecting dogs and cats. While a diagnosis can initially feel overwhelming for pet owners, appropriate treatment, regular monitoring, and strong partnerships between clinics and clients can help many diabetic pets enjoy a good quality of life.

Successful diabetes management relies on an individualized treatment plan, client education, accurate medication administration, appropriate nutrition, and ongoing communication. Veterinary guidelines emphasize the importance of helping owners understand and carry out their pet's care plan at home.

Helping Clients Build Confidence

Insulin injections can initially feel intimidating for owners. Demonstrating proper injection technique, discussing storage and handling of insulin, and ensuring clients understand the importance of consistent dosing schedules can improve confidence and treatment compliance. Current AAHA recommendations also reinforce the importance of matching the correct syringe to the insulin concentration (U-40 or U-100) and selecting syringe sizes that allow for accurate measurement of small doses.

Even small improvements in ease of administration can make a meaningful difference for owners who are committing to twice-daily injections over many years.

Accuracy Matters

Insulin is a medication where precision counts. Small dosing errors may result in poor glycemic control or hypoglycemia, making accurate dose measurement essential.

For patients receiving lower insulin doses, veterinary guidelines recommend insulin syringes that improve dose visualization, helping owners measure doses more confidently and consistently.

Elimedical Partnering with Veterinary Teams

As a Canadian supplier to veterinary practices, Elimedical is committed to supporting clinics with dependable diabetic care supplies designed to help make insulin administration as straightforward as possible.

Elimedical line of Insulin Syringes feature:



U-40 and
U-100 options



0.3 mL, 0.5 mL,
and 1.0 mL
capacities



Fine-gauge
needles designed
for patient comfort



Clear, easy-to-read
graduations to
support accurate
dosing



Reliable quality for
everyday clinical
and at-home use



12 Sizes to Choose From:

Choose the size that is right for you and your patient.

Elimedical Code	Size	Units
03V3108	0.3cc U40, 31G x 5/16"	100pcs/box
03V3108B	0.3cc U100, 31G x 5/16"	100pcs/box
03V2913B	0.3cc U100, 29G x 1/2"	100pcs/box
05V2713VB	0.5cc U100, 27G x 1/2"	100pcs/box
05V2913B	0.5cc U40, 29G x 1/2"	100pcs/box
05V2913VB	0.5cc U100, 29G x 1/2"	100pcs/box
10LDS2713B	1cc U100, 27G x 1/2"	100pcs/box
10LDS2913B	1cc U100, 29G x 1/2"	100pcs/box
10V2813BB	1cc U40, 28G x 1/2"	100pcs/box
10V2813B	1cc U40, 28G x 1/2"	30 pcs/box
40U03M2913	0.3cc U40, 29G x 1/2"	100pcs/box
40U10M2913	1cc U40, 29G x 1/2"	100pcs/box

Whether supporting a newly diagnosed patient or a long-term diabetic, having the right supplies helps reinforce the confidence that both veterinary teams and pet owners need.

This Pet Diabetes Month, Elimedical is proud to support veterinary professionals in helping diabetic patients receive the precise, consistent care they deserve.

